




LED PDT SKIN SKIN THERAPY MACHINE

Germany import LED diode • Seven different colors • Even Energy Distribution



- Wrinkles
- Scarring
- Pigmentation
- Tightening
- Acne
- Skin Rejuvenation
- Rosacea



KLINIKK ANGEL

Aesthetics & Beauty Equipment



PRINCIPLE OF TREATMENT

LED Therapy is the therapeutic use of light to modulate cellular function to improve and quicken wound healing, treat acne, rejuvenate the appearance of skin, stimulate hair growth, improve local circulation, perform 5-ALA photodynamic therapy (PDT), and relieve pain and stiffness in muscles and joints.

The light is placed over the skin with no contact and is illuminated for a period of 15-30 minutes. This allows absorption of photons (light particles) into the target cellular components, resulting in the production of new cellular energy. The response of the inflammatory-stage cells is enhanced, and small proteins are released which support the growth, survival and differentiation of new cells.

7 DIFFERENT COLORS



Purple Light

400-420nm

A combination of red and blue light creates purple light, which helps eliminate toxins, reduces inflammation and improves healing and circulation of blood.

Blue Light

430-450nm

Reduces sebaceous activity and acne vulgaris. Works to clear congestion and kills acne-causing bacteria.

CyanLight

500-520nm

Calms and soothes irritated skin. Supports in reducing the size of swollen capillaries.

Yellow Light

570-600nm

Stimulates collagen production and assists wound repair to naturally heal and regenerate skin from within. Beneficial for rosacea and post-laser treatments.

Orange Light

590-630nm

Promotes a brighter glow for more radiant skin. Adds vitality to dull, lifeless complexions.

White Light


(mixed light)400-830nm

Accelerate active tissue metabolism, reduce fine lines and sagging skin. Rapid cell turn over

Red Light

630-830nm

Increases blood circulation and lymphatic flow to encourage cellular renewal and promote a youthful complexion.




ADVANTAGES



USER FRIENDLY OPERATION INTERFACE



Flexible arm can adjust to your treatment area.



Lifting range: 25CM




Mechanical spring structure,
long life and more stable



180 degree free rotation

Adjustable head for large coverage or focused intensity.



Even Energy Distribution



Optimised Flat Distribution (OFT) allows even and consistent energy distribution throughout the entire target layer and area resulting in consistent results, every time.



Germany imported
high-power LED,
100mw/pc.

LED LIGHT THERAPY BENEFITS INCLUDE

- + Repair and rejuvenate skin
- + Improve skin tone and texture
- + Decreases deeper lines and wrinkles
- + Increases skin's natural hydration
- + Increases blood flow to tissue
- + Decreases the risk of bruising
- + Stimulates ATD - cellular energy
- + Wound Healing
- + Stimulates collagen production
- + Treat acne, dermatitis, rosacea and psoriasis
- + Skin tightening
- + Stimulates serotonin to the brain
- + Relief of muscle and joint pain/pain management
- + Decreases recovery time



SPECIFICATION

Model	PDT LIGHT
Diode Quantity	1500 / 3000
Power	1000W
Diode Power	100mw/diode
Wavelength	400-830nm
LED Colors	Red Blue Yellow Orange Cyan Purple White
Screen	8inch touch screen
Working Mode	Pulse mode/ Continue Mode



KLINIKK ANGEL

Aesthetics & Beauty Equipment

Melhusvegen 505, 7224 Melhus, Thoragården Building, 2 Floor
Mob: +47 979 789 41 Tel : +47 416 963 95
Email: angel@klinikkangel.no

